## QualiSense ADHESIVE S2

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## QualiSense ADHESIVE S2

### Self-etching light cure bonding system

#### Instructions for Use

**ADHESIVE S2** is a simple to use self-etching light curing bonding system. It consists of a self-etching primer **Part A** and a bonding agent **Part B**, that are applied one after the other.

**ADHESIVE S2** designed for bonding of light cure composites to enamel and dentine, silicate ceramics, zirconia and non-precious metals. **ADHESIVE S2** is compatible with all visible light cure composite restorative materials

#### Indications

Adhesive for:

- Composite based restorations with light cure composites

#### Side effects

In singular cases, **ADHESIVE S2** may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material.

#### Incompatibility with Other Materials

Do **not** use in combination with substances containing eugenol because eugenol inhibits the polymerization of **ADHESIVE S2**. Neither store the material in proximity of eugenol containing products, nor let the material allow coming into contact with materials containing eugenol.

#### Application

##### 1. Isolation

Rubber dam is the recommended method of isolation.

##### 2. Cavity Preparation

Clean the tooth with flour of pumice and water prior to preparation. Prepare the cavity with minimal tooth reduction. Margins should have a slight (0.5 - 1.0 mm) bevel placed in the enamel to increase the surface area for greater bond strength.

##### 3. Pulp Protection

For pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material (e.g. **CH LINER**).

##### 4. Application of ADHESIVE S2

###### 4.1. Application of Part A

Apply **Part A** using a brush onto the slightly wet enamel and dentin surfaces (as they are usual in clinical procedures) for about 10 seconds with agitation.

The material should build a homogeneous layer. Air thin gently to remove the volatile components and to disperse the adhesive.

###### 4.2. Application of Part B

Apply **Part B** using a brush for about 10 seconds with agitation. Air thin to remove the volatile components and to disperse the bonding to an even layer until the dentine surface appears slightly dull (the enamel surface remains glossy). Then light cure with a suitable dental halogen light unit or an LED (wavelength 400–500 nm, light intensity min. 1000 mW/cm<sup>2</sup>) for 20 seconds before placement of a light cure composite.

#### Note:

Application of **Part A** and **Part B** can be done with one brush.

##### 5. Restorative Placement

Apply the restorative material according to the instructions of the manufacturer.

Best results are obtained with application of a thin layer of a light cure flowable composite followed by the application of a moldable composite. Light cure each composite layer separately according to the corresponding user instructions.

#### Warnings

- Unpolymerized material may have an irritating effect and may lead to a sensitizing reaction against methacrylates
- Avoid contact with skin, mucous membrane and eyes
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

#### Storage

Do **not** store above 25 °C (77 °F)! Avoid storage in direct sunlight. Do **not** use after expiration date.

#### Warranty

DS Dental Supply GmbH warrants this product will be free from defects

in material and manufacture. DS Dental Supply GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and DS Dental Supply GmbH's sole obligation shall be repair or replacement of the DS Dental Supply GmbH product.

#### Limitation of Liability

Except where prohibited by law, DS Dental Supply GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

#### Keep away from children!

#### For dental use only!