QualiSense LC DESENSITIZER

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QualiSense LC DESENSITIZER

Light Curing One Bottle Desensitizing & Adhesive Agent

Product description

LC DESENSITIZER is a light curing dentine adhesive agent, which serves all at once as a desensitizer. It is designed for strong bonding of **composites, compomers** and metals to **enamel and dentine**. The application as a desensitizer in cervical regions and free dentine areas during and after dental treatments anodynes all postoperative pains. LC DESENSITIZER does not tan and irritate the tissue and gingiva.

The strong adhesion of LC DESENSITIZER to etched enamel or dentine proceeds on principles similar to that occurring with glass ionomer cements. Good, long lasting adhesive strength and good biocompatibility are attained by methacrylate grafted polycarboxylic acids.

The desensitizing LC DESENSITIZER has good wetting properties. It is ethanol-based and hydrophilic. Therefore it can be used on slightly moist dentine surfaces (wet-bonding-technique).

LC DESENSITIZER contains 1000 ppm (0,1 %) fluoride.

Indications/Intended use

- Adhesive/Desensitizer for composite based restorations with light cure composites
- Application after bleaching procedures
- Treatment of hypersensitive cervical regions

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the material. Neither store the material in proximity of eugenol containing products, nor let the material allow coming into contact with materials containing eugenol.

1. Application before Restorations

1.1. Isolation

Rubber dam is the recommended method of isolation.

1.2. Cavity Preparation

Clean the tooth with flour of pumice and water prior to preparation. Prepare the cavity with minimal tooth reduction. Margins should have a slight (0,5 - 1,0 mm) bevel placed in the enamel to increase the surface area for greater bond strength.

1.3. Pulp Protection

For pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material (e. g. CH LINER).

1.4. Enamel and Dentine Conditioning

Recommended is the total etch technique. Apply ETCHING GEL onto the enamel and dentine surfaces beginning with the enamel bevels. Condition the enamel for at least 15 seconds and the dentine for 15 seconds. Deciduous teeth are etched correspondingly longer. Rinse for 20 seconds with water. Dry it in a water-free and oil-free airstream, but do not desiccate. A slightly wet dentine surface is important for the function of LC DESENSITIZER. The etched enamel bevel should have a chalky white appearance.

Etching Precaution

It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.

1.5. Application of LC DESENSITIZER

Apply LC DESENSITIZER generously with a brush onto the enamel and dentine surfaces for 30 seconds with agitation. The material should build a homogeneous layer. Remove excess material carefully. Dry cautiously with oil free air for about 15 seconds to remove all volatile components and to disperse LC DESENSITIZER to an even layer. Do not desiccate the dentine.

LC DESENSITIZER will not self-care. Cure the LC DESENSITIZER coating by exposing its entire area to a dental halogen light unit for 20 seconds before application of a second layer of LC DESENSITIZER.

Notes:

Do not rinse off the LC DESENSITIZER! If not used immediately, place

dispensed mixture in subdued light to prevent premature polymerization by incident light.

1.6. Application of a second layer of LC DESENSITIZER

Apply again LC DESENSITIZER generously with a brush onto the adhesive surfaces as described under 1.5. before placement of a light cure composite or an amalgam restoration.

Note:

Because light curing under amalgam restorations is impossible, and the primer does not selfcure it is essential to cure the coating by exposing its entire area to a dental halogen light unit for 20 seconds before application of the amalgam filling.

1.7. Restorative placement

Refer to manufacturer instructions of restorative materials.

2. Application after bleaching procedures, scaling and in other cases of hypersensitive areas

2.1. Preparation

Clean the sensitive areas and rinse thoroughly. Dry the sensitive area with a cotton pellet or roll to remove excess moisture. Do not overdry!

2.2. Application of LC DESENSITIZER

Apply LC DESENSITIZER generously onto the enamel and dentine surfaces for 30 seconds with agitation. The material should build a homogeneous layer. Remove excess material carefully. Dry cautiously with oil free air for about 15 seconds to remove all volatile components and to disperse the adhesive to an even layer. Do not desiccate the dentine.

Cure the LC DESENSITIZER coating by exposing its entire area to a dental halogen light unit for 20 seconds before application of a second layer of LC DESENSITIZER. Apply again LC DESENSITIZER as described above and light cure for 20 seconds.

Storage

Do not store above 25 °C (77 °F)!

Avoid storage in direct sunlight.

Do not use after expiration date.

Warnings

- Unpolymerized material may have an irritating effect and may lead to a sensitizing reaction against methacrylates.
- Avoid contact with skin, mucous membrane and eyes.
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Warranty

DS Dental Supply GmbH warrants this product will be free from defects in material and manufacture. DS Dental Supply GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and DS Dental Supply GmbH's sole obligation shall be repair or replacement of the DS Dental Supply GmbH product.

Limitation of Liability

Except where prohibited by law, DS Dental Supply GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.